



Valor Peptides

SAMPLE NAME:
SAMPLE CODE:
RECEIVED DATE:
PUBLISHED DATE:

GLOW
XH9P-W7NQ
09/14/2026
09/18/2026



DECLARED IDENTITY GHK-Cu, BPC-157, TB500 (17-23 Fragment)
MATRIX TYPE: PEPTIDE BLEND
SAMPLE NAME: GLOW
SAMPLE SIZE: 70.0 mg
LOT CODE: ARC-GL70-008




72.09 mg
Total Quantity



99.88%
Total Blend Purity

Sample	Test	Specification	Result	Status
Peptide Blend	IDENTITY	GHK-Cu, BPC-157, TB500 (17-23 Fragment)	GHK-Cu, BPC-157, TB500 (17-23 Fragment)	CONFORMS
Peptide Blend	QUANTITY	MEASURE	72.09 mg	MEASURED
Peptide Blend	PURITY	> 98.0%	99.88%	CONFORMS

HPLC Chromatogram Report



RESULTS REVIEWED BY:

Anthony Burke PhD
Lab Manager
09/18/2026

RESULTS CERTIFIED BY:

Dr. Levi Fried MD
Lab Director
09/18/2026

ACCUMARKLABS.COM
XH9P-W7NQ
09/18/2026





Valor Peptides

SAMPLE NAME:
SAMPLE CODE:
RECEIVED DATE:
PUBLISHED DATE:

GLOW
XH9P-W7NQ
09/14/2026
09/18/2026

Analyte	Test	Specification	Result	Status
GHK-Cu	IDENTITY	GHK-Cu	GHK-Cu	CONFORMS
GHK-Cu	QUANTITY	MEASURE	52.243 mg	MEASURED
GHK-Cu	PURITY	> 98.0%	99.933%	CONFORMS
BPC-157	IDENTITY	BPC-157	BPC-157	CONFORMS
BPC-157	QUANTITY	MEASURE	9.453 mg	MEASURED
BPC-157	PURITY	> 98.0%	99.886%	CONFORMS
TB500 (17-23 Fragment)	IDENTITY	TB500 (17-23 Fragment)	TB500 (17-23 Fragment)	CONFORMS
TB500 (17-23 Fragment)	QUANTITY	MEASURE	10.403 mg	MEASURED
TB500 (17-23 Fragment)	PURITY	> 98.0%	99.587%	CONFORMS

RESULTS CERTIFIED BY:

RESULTS CERTIFIED BY:

Anthony Burke
Anthony Burke PhD
Lab Manager
09/18/2026

Levi Fried
Dr. Levi Fried MD
Lab Director
09/18/2026

ACCUMARKLABS.COM
XH9P-W7NQ
09/18/2026





Valor Peptides

SAMPLE NAME:
SAMPLE CODE:
RECEIVED DATE:
PUBLISHED DATE:

GLOW
XH9P-W7NQ
09/14/2026
09/18/2026

Endotoxin USP85 LAL

Test	Specification (EU/mL)	Result (EU/mL)	LOQ (EU/mL)	Verdict
Endotoxin USP85 LAL	≤5	< LOQ	2.5	Conforms

Method: Endotoxin testing performed using Limulus Amebocyte Lysate assay in accordance with USP <85> under validated laboratory conditions.

Rapid Sterility Screening (PCR)

Test	Specification	Result	Unit	Verdict
Rapid Sterility Screening (PCR)	Not Detected	Not Detected		Conforms

Method: Rapid sterility screening performed by real-time polymerase chain reaction (qPCR) targeting conserved microbial DNA sequences under validated laboratory conditions.

RESULTS CERTIFIED BY:


Anthony Burke PhD
Lab Manager
09/18/2026

RESULTS CERTIFIED BY:


Dr. Levi Fried MD
Lab Director
09/18/2026

ACCUMARKLABS.COM
XH9P-W7NQ
09/18/2026

