



Valor Peptides

SAMPLE NAME:
SAMPLE CODE:
RECEIVED DATE:
PUBLISHED DATE:

BPC-157 + TB500 (17-23
Fragment)
MLVW-YB2R
09/15/2026
09/22/2026



DECLARED IDENTITY BPC-157, TB500 (17-23 Fragment)

MATRIX TYPE: PEPTIDE BLEND

SAMPLE NAME: BPC-157 + TB500 (17-23 Fragment)

SAMPLE SIZE: 60.0 mg

LOT CODE: ARC-WV60-001




60.2 mg
Total Quantity



99.91%
Total Blend Purity

Sample	Test	Specification	Result	Status
Peptide Blend	IDENTITY	BPC-157, TB500 (17-23 Fragment)	BPC-157, TB500 (17-23 Fragment)	CONFORMS
Peptide Blend	QUANTITY	MEASURE	60.2 mg	MEASURED
Peptide Blend	PURITY	> 98.0%	99.91%	CONFORMS

HPLC Chromatogram Report



RESULTS REVIEWED BY:

Anthony Burke PhD
Lab Manager
09/22/2026

RESULTS CERTIFIED BY:

Dr. Levi Fried MD
Lab Director
09/22/2026

ACCUMARKLABS.COM
MLVW-YB2R
09/22/2026





Valor Peptides

SAMPLE NAME:
SAMPLE CODE:
RECEIVED DATE:
PUBLISHED DATE:

BPC-157 + TB500 (17-23
Fragment)
MLVW-YB2R
09/15/2026
09/22/2026

Analyte	Test	Specification	Result	Status
BPC-157	IDENTITY	BPC-157	BPC-157	CONFORMS
BPC-157	QUANTITY	MEASURE	29.463 mg	MEASURED
BPC-157	PURITY	> 98.0%	99.848%	CONFORMS
TB500 (17-23 Fragment)	IDENTITY	TB500 (17-23 Fragment)	TB500 (17-23 Fragment)	CONFORMS
TB500 (17-23 Fragment)	QUANTITY	MEASURE	30.743 mg	MEASURED
TB500 (17-23 Fragment)	PURITY	> 98.0%	99.968%	CONFORMS

RESULTS CERTIFIED BY:

RESULTS CERTIFIED BY:


Anthony Burke PhD
Lab Manager
09/22/2026


Dr. Levi Fried MD
Lab Director
09/22/2026

ACCUMARKLABS.COM
MLVW-YB2R
09/22/2026





Valor Peptides

SAMPLE NAME:
SAMPLE CODE:
RECEIVED DATE:
PUBLISHED DATE:

BPC-157 + TB500 (17-23
Fragment)
MLVW-YB2R
09/15/2026
09/22/2026

Endotoxin USP85 LAL

Test	Specification (EU/mL)	Result (EU/mL)	LOQ (EU/mL)	Verdict
Endotoxin USP85 LAL	≤5	< LOQ	2.5	Conforms

Method: Endotoxin testing performed using Limulus Amebocyte Lysate assay in accordance with USP <85> under validated laboratory conditions.

Rapid Sterility Screening (PCR)

Test	Specification	Result	Unit	Verdict
Rapid Sterility Screening (PCR)	Not Detected	Not Detected		Conforms

Method: Rapid sterility screening performed by real-time polymerase chain reaction (qPCR) targeting conserved microbial DNA sequences under validated laboratory conditions.

RESULTS CERTIFIED BY:

RESULTS CERTIFIED BY:


Anthony Burke PhD
Lab Manager
09/22/2026


Dr. Levi Fried MD
Lab Director
09/22/2026

ACCUMARKLABS.COM
MLVW-YB2R
09/22/2026

